



Account Name: \_\_\_\_\_ Account Number: \_\_\_\_\_

Clinic Name: \_\_\_\_\_

Address: \_\_\_\_\_

**Fall Protection  
Brace**

City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_

Prescription  
Order Form

Phone: \_\_\_\_\_ Email: \_\_\_\_\_

Patient Name: \_\_\_\_\_ Age: \_\_\_\_\_ Ht: \_\_\_\_\_ Wt: \_\_\_\_\_

Shoe Type/Size: \_\_\_\_\_ Shoes Enclosed: Yes:  No:  Male:  Female:

Diagnosis, Treatment Goals and Special Instructions \_\_\_\_\_  2-day in-house Rush \$65  3-day in-house Rush \$50

Describe patient condition:

- Left
- Right
- Bilateral

**Treatment Goals:**

- Improve mobility
- Improve lower extremity stability
- Decrease pain \_\_\_\_\_
- Facilitate muscular coordination and gait stability
- Reduce postural sway and increase ankle stability

Color Options

- Black
- Diabetic

Fall Protection Brace  9" (std)  8"  7"  6" (The brace cannot be made higher than the cast)

**Fall Protection Brace Modifications:**

**Top cover**

Polyester cloth covering  
 With soft interface

**Width:**

Normal (std)  
 Narrow  
 Wide

**Length:**

Mets (std)  
 Sulcus  
 Toes

**Shell Options: Polypropylene**

3mm (std)  
 4mm More Rigid  
 5mm Extra Rigid

**Closure:**

Velcro

**Rearfoot**

Intrinsic Post \_\_\_\_\_ °  
 No Post  
 Stable Heel

**Medial Heel Skive:**

2mm  
 4mm  
 6mm

**Calcaneus Alignment**

\_\_\_\_\_ ° Inverted  
 \_\_\_\_\_ ° Everted

**Leg Alignment to the Floor**

\_\_\_\_\_ ° Varum  
 \_\_\_\_\_ ° Valgum

**Footplate Accommodation:**

Navicular  
 Fascia Band  
 Styloid 5th Met  
 Other: \_\_\_\_\_

**Footplate Modifications:**

Medial Flange  
 Lateral Flange

**Other Modifications:**

L  R Heel Lift  1mm  2mm  3mm  
 Medial Flip  
 Arch Pad  1/16"  1/8"  1/4"

**Forefoot Accommodations:**

Balance Pad  L \_\_\_\_\_  R \_\_\_\_\_  
 IPJ Accom  L \_\_\_\_\_  R \_\_\_\_\_  
 Toe Filler  L \_\_\_\_\_  R \_\_\_\_\_  
 Toe Crest  L \_\_\_\_\_  R \_\_\_\_\_  
 Met Pad  L \_\_\_\_\_  R \_\_\_\_\_  
 Met Bar  L \_\_\_\_\_  R \_\_\_\_\_

**Heel Accommodations:**

Heel Spur Donut  1/16"  1/8"  1/4"  
 Heel Spur Horseshoe  1/16"  1/8"  1/4"  
 Heel Pad  1/16"  1/8"  1/4"

Morton's Extension  Enduro  Corex  Plastazote  Shell  Reverse

Check here for a courtesy (Medium) Midleg STS casting sock returned with this order



6299 Guion Rd., Indianapolis, IN 46268 p. 800-444-3632 f. 800-233-2280  
 www.alliedosilabs.com clientservice@aolabs.com