

Caring For Your Child's Custom Orthotics

- ⦿ Wipe clean with a damp cloth. **Avoid soaking or machine washing.**
- ⦿ Check monthly for wear, damage, or fit issues. We recommend purchasing the **Assured Orthotic Replacement Coverage** to protect your against the unforeseen replacement costs due to outgrowth, accidental destruction, and lost or stolen orthotics.
- ⦿ If your child gets new shoes, make sure to transfer your orthotics over; **consult your podiatrist for proper shoe fitting.**
- ⦿ **Do not attempt to trim or alter your child's orthotics.** Adjustments are available and should only be done by your podiatrist or Allied OSI Labs.

Follow-Up Schedule / Checklist

- ☐ **One Month** From Orthotic Fitting
- ☐ **6-Month** Growth & Gait Assessment
- ☐ **1-Year** Rescan

Contact Your Podiatrist Before Scheduled Follow-Ups If:

- ⦿ Your child displays or complains about discomfort, or if there is any swelling and/or redness.
- ⦿ You notice any new instability or increased tripping.
- ⦿ Your child has a sudden growth spurt.

ADD PRACTICE
CONTACT
INFORMATION
HERE