

CUSTOM ORTHOTICS

Accommodative Line

Accommodative

Providing Correction, Control,
& Support With Flexibility



Diabetic Soft

Shell: 5/8" white EVA

Heel Cup Depth: 12mm


Forefoot Post: Intrinsic

Rearfoot Post: Extrinsic EVA

Length: Full

Top Cover: 1/8" plastazote; 1/8"
Poron® laminated

Bottom Cover: 1/16" EVA



Diabetic Medium

Shell: Thermocork

Heel Cup Depth: 12mm

Forefoot Post: Intrinsic

Rearfoot Post: Extrinsic thermocork

Length: Full

Top Cover: 1/8" plastazote; 1/8"
Poron® laminated

Bottom Cover: 1/16" EVA



Diabetic Firm

Shell: 2mm polyethylene

Heel Cup Depth: 16mm

Forefoot Post: Intrinsic

Rearfoot Post: Extrinsic crepe

Length: Full

Top Cover: 1/8" plastazote; 1/8"
Poron® laminated

Bottom Cover: 1/16" EVA

Features: Medial flange, lateral clip, and
deep heel cup



Classic Leather

Shell: Leather, THK, and corax

Heel Cup Depth: 12mm

Forefoot Post: None

Rearfoot Post: Extrinsic corax

Length: Sulcus

Padding: 1/16" plastazote

Top Cover: Leather

Bottom Cover: Ultrasuede®